

## A. Reagents Preparation

### Transfer Buffer:

3.03 g Tris  
14.4 g Glycine  
1 g SDS  
200 ml Methanol  
500 ml MilliQ  
Adjust **pH 8.3-8.5**  
Complete to 1L MilliQ

Proteomics & Metabolomics Unit  
[proteomics.lab@57357.org](mailto:proteomics.lab@57357.org)

Children's Cancer Hospital 57357  
Cairo, Egypt

<https://www.57357.org/proteomics-unit/>

### 10x TBS:

24 g Tris base  
88 g NaCl  
700 ml MilliQ  
Adjust **pH 7.6**  
Complete to 1L MilliQ

### 1x TBST: (0.1% Tween)

100 ml 10x TBS  
900 ml MilliQ  
1 ml Tween 20

### Blocking Solution (Milk):

5% non-fat milk in TBST

### Blocking Solution (BSA):

5% BSA in TBST

### Medium Stripping Solution:


15 g Glycine  
1 g SDS  
10 ml Tween 20

### High Stripping solution:

20 ml 10% SDS  
12.5 ml Tris HCl pH 6.8  
67.5 ml MilliQ  
0.8 ml BME

## B. Transfer:

1. Soak 4 whatman filter papers in cold transfer buffer for 20 min
2. Soak membrane in methanol in 30 sec. Discard. Then soak in cold transfer buffer for 30 min while shaking. (You could change the transfer buffer one time to make sure methanol was washed)
3. Soak the gel (10%) in running buffer after removing it from the cassette.
4. Assemble the transfer sandwich in the semi-dry trans-blotter as follows:



5. Close the trans-blotter and run at 25V constant, 1.0 A limit, 10 minutes, using Thermo/BioRad blotter\*.

\*Run conditions depend on the used antibody.

## C. Immunoblotting part I:

6. Stain the gel in coomassie stain for ~ 20 min. Then remove stain and add destain solution
7. Wash the membrane for 2 min in TBST while shaking.
8. Prepare either 3-5% non-fat dry milk in TBST solution or 3-5% BSA in TBST solution.
9. Remove the TBST solution from the membrane.
10. For **Blocking**: Incubate membrane in 3-5% non-fat dry milk or 3-5% BSA for 1 hour at room temperature with shaking.

Proteomics & Metabolomics Unit  
[proteomics.lab@57357.org](mailto:proteomics.lab@57357.org)

Children's Cancer Hospital 57357  
Cairo, Egypt

<https://www.57357.org/proteomics-unit/>

11. Prepare 3 ml of primary antibody dilution in either 3-5% non-fat dry milk or 3-5% BSA as indicated in the antibody's data sheet.
12. Remove the blocking solution from the membrane and incubate the membrane in primary antibody solution overnight at 4°C with shaking.

**D. Immunoblotting part II:**

1. Remove primary antibody solution.
2. Wash membrane with 1X TBST 3 times, 15 minutes each.
3. Prepare 3 ml secondary antibody dilution as indicated in its data sheet in 3-5% non-fat dry milk solution or 3-5% BSA in TBST.
4. Remove TBST from the membrane and incubate the membrane in secondary antibody solution at room temperature for 1 hour with shaking.

**E. Detection:**

1. Remove secondary antibody solution.
2. Wash membrane with 1X TBST 3-4 times, 5 minutes each.
3. Prepare ECL solution by mixing both components at a ratio of 1:1 prior to detection (~2 ml per membrane).
4. Drain TBS from the membrane and incubate the membrane in ECL solution for 5 minutes at room temperature.
5. Drain excess ECL and visualize bands using chemidoc detection system.

Proteomics & Metabolomics Unit  
[proteomics.lab@57357.org](mailto:proteomics.lab@57357.org)

Children's Cancer Hospital 57357  
Cairo, Egypt

<https://www.57357.org/proteomics-unit/>